

ALLAMERICAN DENTAL LAB

Implant Brand	Platform Type	Titanium	Zirconia	Hybrid	Screw Torque (N-CM)
Astra Tech	OsseoSpeed 3.5/4.0	✓		✓	25
	OsseoSpeed 4.5/5.0	✓		✓	25
	Fixture MT 3.5	✓			25
	Fixture MT 4.0	✓			25
	Fixture MT 4.5/5.0	✓			25
BioHorizons®	External 3.5	✓	✓		30
	External 4.0	✓	✓		30
	External 5.0	✓	✓		30
	External 6.0	✓	✓		30
	Internal 3.0	✓		✓	30
	Internal 3.5	✓		✓	30
	Internal 4.5	✓		✓	30
	Internal 5.7	✓		✓	30
Biomet 3i™	Certain 3.4	✓		✓	20
	Certain 4.1	✓		✓	20
	Certain 5.0	✓		✓	20
	Certain 6.0	✓		✓	20
	External 3.4	✓	✓		32-35
	External 4.1	✓	✓		32-35
	External 5.0	✓	✓		32-35
	External 6.0	✓	✓		32-35
Camlog®	Φ3.3	✓		✓	20
	Φ3.8	✓		✓	20
	Φ4.3	✓		✓	20
	Φ5.0	✓		✓	20
	Φ6.0	✓		✓	20
Conelog®	Φ3.3	✓		✓	20
	Φ3.8	✓		✓	20
	Φ4.3	✓		✓	20
	Φ5.0	✓		✓	20
Friadent XiVE®	Frialit 3.0	✓			24
	Frialit 3.4	✓		✓	24
	Frialit 3.8	✓		✓	24
	Frialit 4.5	✓		✓	24
	Frialit 5.5	✓		✓	24
	Frialit 6.5	✓		✓	24
Intra-Lock®	SQ (Straight 3.3/Taper 3.75)	✓			30
	Standard (Straight 4.0/Conic 4.3)	✓		✓	30
	Wide (Straight 4.75/Conic 4.3, 5.0)	✓		✓	30
Implant One Platform	Series 100	✓		✓	20
	Series 200	✓		✓	20



✓ = Available
 For 1-piece abutments, either brand screws or All American screws may be used. Implant brand screws are not provided by All American

All registered marks are owned by the respective manufacturers and are not the property of All American

ALLAMERICAN DENTAL LAB

Implant Brand	Platform Type	Titanium	Zirconia	Hybrid	Screw Torque (N-CM)
LifeCore	External Hex SD	✓	✓		30
	External Hex RD	✓	✓		30
	External Hex WD	✓	✓		30
	Internal Hex SDI	✓		✓	30
	Internal Hex RDI 4.5/4.75	✓		✓	30
Megagen <small>*Non Engaging Options Available</small>	AnyRidge Internal *			✓	30
	AnyOne Internal *			✓	30
	EZ Plus Internal Small			✓	30
	EZ Plus Internal Reg.			✓	30
	ExFeel Internal			✓	30
	Rescue Internal			✓	30
Neoss®	3.25	✓		✓	32
	3.5-5.5	✓		✓	32
Nobel Biocare®	Branemark NP	✓	✓		35
	Branemark RP	✓	✓		35
	Branemark WP	✓	✓		35
	Replace NP	✓		✓	35
	Replace RP	✓		✓	35
	Replace WP	✓		✓	35
	Replace 6.0	✓		✓	35
	Active NP	✓		✓	35
	Active RP	✓		✓	35
Osstem®	GS Mini (Φ3.5)	✓		✓	30
	GS Standard (Φ4.0/Φ4.5/Φ5.0)	✓		✓	30
	SS Regular (Φ4.8)	✓		✓	30
	SS Wide (Φ6.0)	✓		✓	30
	US Mini (Φ3.5)	✓			30
	US Regular (Φ4.1)	✓			30
	US Wide (Φ5.1)	✓			30
Straumann®	synOcta® ITI NN	✓	✓		35
	synOcta® ITI RN	✓		✓	35
	synOcta® ITI WN	✓		✓	35
	Bonelevel NC (Φ3.3)	✓		✓	35
	Bonelevel RC (Φ4.1/Φ4.8)	✓		✓	35
Zimmer®	Tapered Screw-vent 3.5	✓		✓	30
	Tapered Screw-vent 4.5	✓		✓	30
	Tapered Screw-vent 5.7	✓		✓	30
	Spline 3.25	✓	✓		28.2
	Spline 3.75/4.0	✓	✓		28.2
	Spline 5.0	✓	✓		28.2

✓ = Available

For 1 -piece abutments, either brand screws or All American screws may be used. Implant brand screws are not provided by All American

All registered marks are owned by the respective manufacturers and are not the property of All American

